

Find differences in protein expression with Ettan DIGE

2-D Fluorescence Difference Gel Electrophoresis



Results you can trust in protein abundance analysis

Ettan™ DIGE system uses 2-D Fluorescence Difference Gel Electrophoresis (2-D DIGE), the benchmark for protein abundance analysis. The world's leading pharmaceutical and academic centers of excellence are increasingly switching to Ettan DIGE because of the unrivaled accuracy it provides.

Standardize your 2-D results


Ettan DIGE uses multiplexing, the simultaneous co-separation of multiple, fluorescently labeled samples, including a pooled internal standard on each gel. Each protein spot has its own internal standard, which ensures that the differences you see in protein abundance are real. This is the only effective way to minimize gel-to-gel variation and significantly increase accuracy and reproducibility.

Efficiency saves time and costs

The level of reproducibility and statistical confidence with Ettan DIGE ensures that you get results you can rely on. You can also be confident that you're not missing important differences in protein abundance. Because Ettan DIGE requires far fewer gels than other methods for measuring protein abundance differences, it also provides significant time and cost savings.

Ettan DIGE delivers these important benefits:

- **Standardization:** uses a pooled internal standard to deliver the lowest possible gel-to-gel variation
- **Accuracy:** ensures that the smallest possible real differences in protein abundance are detected with unparalleled statistical confidence
- **Reproducibility:** leads to consistent and comparable conclusions
- **Cost effectiveness:** achieves dependable results with far fewer gels
- **Efficiency:** increases throughput and significantly reduces analysis time
- **Proven method:** has rapidly become established in leading laboratories worldwide



Ettan DIGE: raising the standard of 2-D electrophoresis

Conclusions you can defend

Accurate results lead to confident conclusions and the right decisions. In proteomics research, you need to be sure you are choosing the right proteins to select, characterize, and quantitate. Ettan DIGE minimizes the danger of missing vital biomarkers, so you can be confident you have selected the right candidates.

Ettan DIGE is a proven method that combines the most statistically accurate protein abundance data with the greatest resolving power available in 2-D electrophoresis. Whether you're interested in researching disease, identifying biomarkers, or understanding tissue differentiation, Ettan DIGE gives you findings you can believe in, which lead to conclusions you can defend.

“Without the benefit of the internal standard, 42 of these proteins would have been overlooked due to variation between normal and tumor samples ... compared with individual DIGE comparisons made within a single 2-D gel.”



Friedman, D. B. *et al.* Proteome analysis of human colon cancer by two-dimensional difference gel electrophoresis and mass spectrometry. *Proteomics* **4**, 793-811 (2004).

The use of a pooled internal standard leads to accurate conclusions

Traditional 2-D electrophoresis

Four different samples run on four different gels

The abundance of this particular protein spot appears to be increasing in samples 3 and 4. Is this increase due to system variation or induced biological change?





Experimental conclusion

Without running a significant number of replicates to average the results, the conclusion would be an increase in abundance in samples 3 and 4.

Ettan DIGE using a pooled internal standard

Four different samples, plus one internal standard, on two different gels



Experimental conclusion

The same internal standard is run on both gels. The increase in abundance of protein in gel B, as shown by the increase in the internal standard, is due to gel-to-gel variation. When the internal standard is normalized between gels A and B, the conclusion is that the abundance of protein in sample 3 has actually decreased.

The differences you see in protein abundance are real

Ettan DIGE is the only system with a pooled internal standard for every spot on every gel. Ettan DIGE uses prelabeling of protein samples with size- and charge-matched, spectrally resolvable CyDye™ DIGE Fluors.

After separation on Ettan IPGphor™ or Multiphor™ II IEF system in the first dimension and Ettan DALTwelve or Ettan DALSix Large Vertical Electrophoresis System in the second dimension, samples are scanned with Typhoon™ Variable Mode Imager or Ettan DIGE Imager.

Differences in protein abundance are then accurately quantitated using DeCyder™ 2-D Differential Analysis Software.

This dedicated 2-D analysis software uses the internal standard to derive data from within gels and then between gels, minimizing gel-to-gel variation. This allows you to achieve detection of < 10% differences between samples with > 95% statistical confidence, within minutes.

“Ettan DIGE with the internal standard delivers reliable results by decreasing gel-to-gel variation more than ever before.”


Dr. Jörgen Östling, AstraZeneca R&D, Mölndal, Sweden

Ettan DIGE performance comparison*

	Matched spots	Mean CV%
DIGE dyes, no internal standard, without DeCyder 2-D software	1300	26%
DIGE dyes, with internal standard, without DeCyder 2-D software	1300	15%
Ettan DIGE, with DeCyder 2-D Differential Analysis Software	1300	3-6%
Classical 2-D with fluorescent stain, without DeCyder 2-D software	1670	21%

* For each method, the same sample was analyzed in replicates of six gels.

Ettan DIGE is the only system with a pooled internal standard for every spot on every gel



Mix labeled extracts

2-D electrophoresis separation

Typhoon Variable Mode Imager or Ettan DIGE Imager

DeCyder 2-D Differential Analysis Software



1. Sample labeling

CyDye DIGE Fluors

Designed specifically for Ettan DIGE, these dyes are size- and charge-matched, as well as spectrally distinct, offering bright and intense colors with narrow excitation and emission bands. This allows co-separation of different CyDye DIGE fluor-labeled samples in the same gel and ensures that all samples are subject to exactly the same first- and second-dimension electrophoresis running conditions. This limits experimental variation and ensures accurate spot matching within gels.

CyDye DIGE Fluor minimal dyes

These CyDye fluors label the amino acid lysine. They are the most commonly used dyes for general 2-D electrophoresis applications and are suitable for all sample types. Each labeling reaction requires 50 µg of protein sample.

CyDye DIGE Fluor Labeling Kit for scarce samples

These CyDye fluors label the amino acid cysteine and are recommended for use when scarce amounts of sample are available (e.g., for analysis of proteins obtained from laser-capture microdissection). Each labeling reaction requires only 5 µg of protein sample.



2. Sample separation: first dimension

Ettan IPGphor 3 and Multiphor II IEF systems

Ettan IPGphor 3 IEF allows for high throughput first dimension electrophoresis with increased speed, accuracy and reproducibility. With Ettan IPGphor Manifold or Ettan IPGphor Manifold Light, up to 12 IPG strips can be run simultaneously which simplifies handling of multiple

samples and saves time. Ettan IPGphor 3 Control Software helps you manage up to four Ettan IPGphor 3 units simultaneously.

To perform first dimension separation on the Multiphor II Electrophoresis System, use the Immobiline™ DryStrip Kit. To increase efficiency and reproducibility even further, the specially designed IPGbox can be

used to rehydrate up to 12 IPG strips independently and at the same time. With the IPGbox, you can load samples during rehydration by including them in the buffer, or by using sample cups.



IPGbox is designed for rehydration of up to 12 Immobiline DryStrip Gels.



Ettan IPGphor 3 system, optimized for speed, reliability and use with fluorescently-labeled protein samples.


3. Sample separation: second dimension

Ettan DALT Electrophoresis System

In just 4-5 hours, Ettan DALTsix Large Vertical System can run up to six large labcast or precast gels (26 x 20 cm) in the second dimension. The gels are compatible with first dimension 18 or 24 cm Immobiline DryStrip Gels (IPG strips) from both Ettan IPGphor and Multiphor II IEF systems. For increased

throughput, use the EttanDALTwelve electrophoresis system which can simultaneously run up to 12 large-format gels. Ettan DALTwelve used with 24 cm IPG strips provides very high resolution and maximizes the number of distinct, quantifiable protein spots per gel.

Choosing a product with the DIGE-approved seal assures that you will obtain optimum results with a product fully compatible with Ettan DIGE technology. From sample preparation to data analysis, there is a DIGE-approved product for every step. Look for the DIGE-approved seal for 2-D results you can trust.



Compact Ettan DALTsix
Large Vertical System



Ettan DALTwelve provides high resolution,
throughput, and batch-to-batch reproducibility.

4. Image acquisition

Typhoon Variable Mode Imager

The Typhoon Variable Mode Imager is a highly sensitive, variable mode gel imager optimized to image CyDye DIGE fluor-labeled proteins.

- Outstanding linearity, quantitative accuracy, and extremely low limits of detection
- Detection of small differences with highest statistical confidence
- Consistent point light illumination
- Simultaneous imaging of two large-format Ettan DALT gels automatically

Ettan DIGE Imager

Ettan DIGE Imager represents leading-edge technology for generating high-quality, low-noise images from your 2-D DIGE experiments. By combining very high resolution with precise motion control, Ettan DIGE Imager produces accurate, multichannel images of your CyTM2-, Cy3-, and Cy5-labeled gel.

- High sensitivity for imaging faint spots
- Software-selectable wavelengths for better sensitivity and reduced background
- Adjustable scan time to optimize sensitivity
- Sealed environment for scanning and protecting wet samples



Gels are scanned between glass plates, preventing drying and shrinkage and allowing further running and rescanning, if required.



Gels are mounted in cassettes to avoid drying.


6. Multivariate analysis in proteomics

DeCyder Extended Data Analysis (EDA) Software


DeCyder EDA is a powerful analysis tool extending the statistical options offered in 2-D DIGE. It enables you to switch seamlessly between EDA statistical results and DeCyder 2-D visualization data, putting your results into biological context by linking to internal and external databases. Together they create a powerful platform of image analysis tools, which offers flexibility, versatility, and adaptability for protein profiling.

DeCyder EDA provides advanced statistical analysis in a simple-to-use format, uncovering patterns in expression data and relationships using multivariate analysis and sophisticated clustering methods.


The software helps you understand regulatory pathways, find proteins with similar expression profiles, or group your samples according to common expression patterns. It also supports the identification of proteins that discriminate between disease stages, tumor types or other sample subtypes, giving your Ettan DIGE results clarity and biological relevance.



Hierarchical cluster analysis



K-means analysis



Principal component analysis (PCA)

Proven in practice

CyDye DIGE Fluor minimal dyes for tracking membrane protein abundance in epilepsy-induced rats

The results of this proteomics study demonstrate that, as part of the Ettan DIGE platform, CyDye DIGE Fluor minimal dyes are applicable to the analysis of membrane proteins and are beneficial for use within preclinical and clinical proteomics.


The study was carried out in collaboration with Dr. S. Hattori *et al.*, Division of Cell Proteomics, Institute of Medical Science, University of Tokyo, Japan.

Some examples of the proteins identified by mass spectrometry

Protein name	Differential expression*
NADH dehydrogenase alpha subunit 10	+ 6.27
Oxygen regulated protein	+ 2.10
ATP-dependent zinc metalloprotease	+ 1.96
Golgi coiled coli protein	+ 1.72
Neural adhesion molecular F3	+ 1.96
Aldehyde dehydrogenase	+ 1.67
Huntingtin interaction protein 1 related protein	+ 4.36
T-cell activation Rho GTPase-activation protein	+ 1.51
Growth factor receptor bound protein	+ 2.34
Androgen-induced prostate proliferative shutoff associated protein As3	+ 2.02
NADH dehydrogenase alpha subunit 10	- 33.02
FERM, RhoGEF and pleckstrin domain protein 2	- 1.67
Heat-stable enterotoxin receptor precursor	- 1.87
ATP synthase beta subunit	- 2.67
F1-ATPase beta subunit	- 2.17
Huntingtin interacting protein related	- 2.63
H ⁺ -transporting two sector ATPase	- 1.87


* Abundance in treated sample vs. control, measured by volume ratio

Analyzing differences in membrane protein expression in rat brain with 2-D DIGE



Protein spots showing a predominantly red or green color demonstrate changes in abundance between normal and kainic acid treatment (used to induce epilepsy in rats). The centrally positioned pair of green and red spots shown above were identified as the same protein, NADH dehydrogenase alpha subcomplex 10-like protein. This protein has undergone some structural changes caused by the kainic acid treatment.

Without 2-D DIGE



Without the 2-D DIGE system, several replicates must be run due to high gel-to-gel variation. This is time-consuming, delivering lower quality data and leading to inaccurate conclusions. Courtesy of Jörgen Östling, AstraZeneca R&D Mölndal, Sweden.

Ordering information

Product	Quantity	Code no.	pH range	Quantity	Code No	
					18 cm	24 cm
Sample labeling			3–5.6 NL	12/pack	17-6003-56	17-6003-57
DIGE Trial Pack, incl. CyDye DIGE Fluor minimal dye labeling kit + 1 month free DeCyder 2-D trial license	2 nmol	28-9373-73	5.3–6.5	12/pack	17-6003-61	17-6003-62
CyDye DIGE Fluor minimal dye labeling kit	2 nmol	28-9345-30	6.2–7.5	12/pack	17-6003-66	17-6003-67
CyDye DIGE Fluor, minimal labeling kit	5 nmol	25-8010-65	7–11 NL	12/pack	17-6003-71	17-6003-72
CyDye DIGE Fluor Labeling Kit for Scarce Samples	1	25-8009-83	3–11 NL	12/pack	17-6003-76	17-6003-77
CyDye DIGE Fluor Labeling Kit for Scarce Samples and Preparative Gel Labeling	1	25-8009-84	¹ IPG strips are also available in 7, 11 and 13 cm strip lengths. For information, contact your local GE Healthcare representative.			
Sample separation, first dimension			Product			
DeStreak Rehydration Solution	5 × 3 ml	17-6003-19	Quantity			
DeStreak Reagent	1 ml	17-6003-18	Code no.			
Ettan IPGphor 3 Isoelectric Focusing Unit	1	11-0033-64	Sample separation, second dimension			
Immobiline DryStrip Kit	1	18-1004-30	Ettan DALTwelve Separation Unit and Power Supply/Control Unit, 115V			
IPGbox	1 IPGbox + IPGbox Kit	28-9334-65	Ettan DALTwelve Separation Unit and Power Supply/Control Unit, 230V			
IPGbox Kit	10 Reswell Trays + IPGbox Insert	28-9334-92	DALTwelve Gel Caster, complete			
IPG Buffer pH 5.5-6.7	1 ml	17-6002-06	Ettan DALTsix Electrophoresis Unit, 115V			
IPG Buffer pH 4-7	1 ml	17-6000-86	Ettan DALTsix Electrophoresis Unit, 220V			
IPG Buffer pH 6-11	1 ml	17-6001-78	EPS 601 Power Supply 18-1130-02			
IPG Buffer pH 7-11 NL	1 ml	17-6004-39	DALTsix Gel Caster			
IPG Buffer pH 3-10 NL	1 ml	17-6000-88	MultiTemp™ III, Thermostatic Circulator, 115 V			
IPG Buffer pH 3-10	1 ml	17-6000-87	MultiTemp III, Thermostatic Circulator, 230 V			
IPG Buffer pH 3-11 NL	1 ml	17-6004-40	Image acquisition			
Ettan IPGphor Manifold, Complete	1	80-6498-38	Ettan DIGE Imager			
Ettan IPGphor Manifold, Light Complete	1	11-0026-88	Typhoon Trio™ with ImageQuant™ TL and PC			
Multiphor II Electrophoresis System	1	18-1018-06	Typhoon Trio+ with ImageQuant TL and PC			

Immobiline DryStrip Gels for IEF¹

pH range	Quantity	Code No	
		18 cm	24 cm
3.5–4.5	12/pack	–	17-6002-38
3–7 NL	12/pack	–	17-6002-43
4–7	12/pack	17-1233-01	17-6002-46
6–9	12/pack	17-6001-88	17-6002-47
6–11	12/pack	17-6001-97	–
3–10	12/pack	17-1234-01	17-6002-44
3–10 NL	12/pack	17-1235-01	17-6002-45

Image analysis		
DeCyder 2-D v7.0 Software Package		28-9435-83
DeCyder 2-D v7.0 Oracle 10gR2, 5 user license	1	28-9435-88
DeCyder 2-D v7.0, 1 user license	1	28-9442-75

Related literature		Code no.
2-D Electrophoresis: Principles and Methods		80-6429-60
Ettan DIGE System User Manual		18-1173-17

For local office contact information, visit,
www.gelifesciences.com/contact

GE Healthcare Bio-Sciences AB
 Björkgatan 30
 751 84 Uppsala
 Sweden

www.gelifesciences.com/dige

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Cy, CyDye, DeCyder, Ettan, ImageQuant, Immobiline, IPGphor, Multiphor, MultiTemp, Typhoon, and Typhoon Trio are trademarks of GE Healthcare companies.

© 2005–2009 General Electric Company –
 All rights reserved. First published Nov. 2005

GE Healthcare UK Limited
 Amersham Place
 Little Chalfont
 Buckinghamshire, HP7 9NA
 UK

GE Healthcare Bio-Sciences Corp.
 800 Centennial Avenue, P.O. Box 1327
 Piscataway, NJ 08855-1327
 USA

GE Healthcare Europe GmbH
 Munzinger Strasse 5, D-79111 Freiburg
 Germany

GE Healthcare Bio-Sciences KK
 Sanken Bldg., 3-25-1, Hyakunincho
 Shinjuku-ku, Tokyo, 169-0073
 Japan

